

MIT OpenCourseWare
<http://ocw.mit.edu>

7.13 Experimental Microbial Genetics

Fall 2008

For information about citing these materials or our Terms of Use, visit: <http://ocw.mit.edu/terms>.

Primer Design Guidelines for PCR Reactions

General types of PCR reactions you may perform:

- Gene amplification (with introduction of restriction sites for downstream cloning)
- Point mutations
- Deletions

General considerations for primer design:

- Primers should be roughly 17-35 bases in length
- T_m (primer melting temperature) ~ 55-70°C
 - most simply calculated as $T_m = [4(G+C) + 2(A+T)]^\circ\text{C}$
 - or use a web-based calculator (e. g. <http://www.idtdna.com/SciTools/SciTools.aspx>)
- Primers should end (3') in a G or C, or CG or GC
 - improves efficiency of priming
- Avoid runs of three or more Cs or Gs at the 3'-ends of primers
 - may result in mis-priming at G or C-rich sequences
- Avoid primer *self*-complementarity (within a single primer such that 2° structures like hairpins form)

```

Self Complementarity
Max complementarity found 6 bp, free energy=-0.30 Kcal/mol

5'-ACGTCCAAGAGGAAGCCCTCTT-3'
*****
A:::C
A::G
G+C
G+C
A+T
G+C
A+T
A+T

5'-ACGTCC-3'
  
```

- Avoid primer *pair* complementarities (primer dimers will form and be synthesised preferentially to any other product)

```

Primer_Primers
Two primers complementarity.
Max complementarity in continuous: 4 bp, free energy=-1.60 Kcal/mol

5'-ACGTCCAAGAGGAAGCG-3'
      ||||
3'-ACTTACGACTTTCGAATGTAA-5'
  
```

- Avoid use of 4 or more consecutive Guanine nucleotides in primer (forms a significantly stable “cruciform” or “guanine tetraplex” 2° structure)

Specific considerations for primer design:

- Remember to maintain plasmid elements required for protein expression:
 - In-frame fusions
 - Stop codons (*reverse primer*)
 - Shine Dalgarno sequence (*forward primer*)
- Always consider restriction sites available in the MCS (multiple cloning site)
- When engineering restriction sites at ends of an amplified PCR product, always include a few additional bases after the restriction site
(http://www.neb.com/nebecomm/tech_reference/restriction_enzymes/cleavage_olignucleotides.asp)

Reference where more details can be found: <http://www.mcb.uct.ac.za/pcroptim.htm>